CoaguChek aPTT Controls

REF	▽	SYSTEM
06882692 190	4 bottles each of lyophilized control plasma Level 1 & 2 <p>8 diluent-filled bulb pipettes (red)</p> <p>1 quality control code chip</p>	CoaguChek® Pro II

English

Intended use

CoaguChek aPTT Controls are used for system checks and quality control testing of activated Partial Thromboplastin Time (aPTT) with the CoaguChek Pro II meter and CoaguChek aPTT Test strips. CoaguChek aPTT Controls are intended for professional use.

Reagents - working solutions

Each bottle of control contains lyophilized anticoagulated human plasma. Each diluent-filled bulb pipette contains calcium chloride and preservatives in water.

Target values and ranges

The code chip provided in the pack contains all relevant lot-specific information to perform a quality control test.

The CoaguChek Pro II meter displays the control range, the reading obtained and an evaluation of the result (inside or outside of the control range). The result is automatically saved as a control reading in the memory of the meter.

The control intervals and limits should be adapted to each laboratory’s individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Precautions and warnings

For in vitro diagnostic use. Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. The testing methods used assays approved by the FDA or cleared in compliance with the European Directive 98/79/EC, Annex II, List A.

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.

Handling

- Allow control bottle and bulb pipette to acclimatize to room temperature at least 15 minutes before reconstitution. Unopened product can be stored at room temperature for the rest of the day.
- Open control bottle.
- Shake all the liquid contained in the bulb pipette out of the capillary into the hollow part of the pipette.
- Use scissors to cut the tip off the bulb pipette. **This should be done 3-5 mm below the welding.** Avoid loss of liquid.
- Add all the liquid in the bulb pipette to the control bottle. Make sure that no liquid is left in the pipette.
- Place the empty pipette into the control vial.
- Dissolve the content of the bottle by gently swirling the bottle a few times. Avoid formation of bubbles.
- During the reconstitution time of 5 minutes periodically swirl the bottle gently a few times. Ensure that the control is completely dissolved.
- Use the control material within 5 minutes of reconstitution.
- Application of sample is performed using the bulb pipette. Aspirate the entire control sample with the pipette before applying a single drop onto the strip.
- Make sure that the drop does not contain air bubbles. If necessary, discard the first drop.

Storage and stability

Store at 2-8 °C.

The lyophilized control plasma is stable up to the stated expiration date.

The reconstituted control material is stable for 5 minutes.

Materials provided

- 4 bottles of lyophilized control plasma Level 1
- 4 bottles of lyophilized control plasma Level 2
- 8 diluent-filled bulb pipettes (red)
- 1 code chip

Materials required (but not provided)

- REF 07210841, CoaguChek Pro II meter (with W-LAN) or
- REF 07237944, CoaguChek Pro II meter (without W-LAN)
- REF 06882382, CoaguChek aPTT Test

Assay

- Prepare the meter and a test strip according to the instructions given in the operator’s manual and in the method sheet for the test strips.
- Insert the quality control code chip into the meter.
- Use the bulb pipette to apply one hanging drop of the control solution to the sample application area of the test strip.
- Use the control solution within 5 minutes after reconstitution.

System checks

The CoaguChek Pro II meter has a number of inbuilt system checks. Please refer to the operator’s manual for details. The CoaguChek aPTT Test strip has an inbuilt quality control function.

Performance data

Reproducibility of CoaguChek aPTT Test was determined using CoaguChek aPTT Controls at 4 external study sites. Controls were measured on 21 days in 2 runs per day and with 2 or 3 lots per study site. Reproducibility figures were calculated by ANOVA (Analysis of Variance).

Reproducibility			
aPTT Control Level	Mean (s)	SD (s)	CV (%)
Level 1	46.6	2.4	5.1
Level 2	66.4	4.4	6.6

For further information, please refer to the appropriate meter manual and the Method Sheets of all necessary components. A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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© 2017, Roche Diagnostics
COAGUCHEK is a trademark of Roche.
All other product names and trademarks are the property of their respective owners.
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim