## Précision

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- Ansell J, Jacobson A, Levy J et al. Guidelines for implementation of patient self-testing and patient self-management of oral anticoagulation. International consensus guidelines prepared by International Self-Monitoring Association for Oral Anticoagulation. Int J Cardiol 2005;99:37-45.
- Menéndez-Jándula B, Souto J et al. Comparing Self-Management of Oral Anticoagulant Therapy with Clinic Management. Ann Intern Med 2005;142:1-10.
- Taborski U, Wittstamm F-J, Bernardo A: Cost-effectiveness of self-managed anticoagulant therapy in Germany. Sem Thromb Hemostas 1999;25:103-109.
- Tripodi A, Breukink-Engbers WGM, van den Besselaar AMHP: Oral Anticoagulant Monitoring by Laboratory or Near-Patient Testing: What a Clinician Should Be Aware Of. Sem Vasc Med 2003;3:243-254.
- Moll S and Ortel TL Monitoring Warfarin Therapy in Patients with Lupus Anticoagulants. Annals of Internal Medicine 1997;127:177-185.

- Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2014.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- Plesch W, Wolf T, Breitenbeck N et al. Results of the performance verification of the CoaguChek XS System. Thromb Res 2008;123:381-389.
- van den Besselaar AMPH, Hoekstra MMCL, Witteveen E, et al. Influence of Blood Collection Systems on the Prothrombin Time and International Sensitivity Index Determined With Human and Rabbit Thromboplastin Reagents. Am J Clin Pathol 127:724-729, 2007.
- van den Besselaar AMPH, Rutten WPF, Witteveen E. Effect of magnesium contamination in evacuated blood collection tubes on the prothrombin time test and ISI calibration using recombinant human thromboplastin and different types of coagulometer.Thrombosis Research 115:239-244, 2005.

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COAGUCHEK and ACCU-CHEK are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2020, Roche Diagnostics

www.coagucheck.com	
	Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim www.Roche.com