CoaguChek XS PT Test PST

[REF]	[▽]	[SYSTEM]
07671679	6	CoaguChek® INRange CoaguChek® XS
07671687	24	
07762798	2x24	

English

Intended use

The CoaguChek XS PT Test PST is an in-vitro assay for the determination of prothrombin time (PT) using the CoaguChek INRange / CoaguChek XS meter. The test uses fresh capillary whole blood and is intended for patient self-testing only.

Test principle

Electrochemical measurement of prothrombin time following activation of blood coagulation with human recombinant tissue factor. Each test strip has a test area containing a prothrombin reagent. When blood is applied, the reagent is dissolved and an electrochemical reaction takes place which is transformed into a clotting time value. The clotting time value is displayed on the meter screen.^{1,2}

These instructions feature two symbols that draw your attention to important information.

⚠ This symbol indicates that the result could possibly be incorrect or that you could be at risk of damaging your health.

ℹ This symbol draws your attention to other important information.

Measuring range and therapeutic range

The CoaguChek INRange/XS meter can display results in the following units: INR (International Normalized Ratio), % Quick and Seconds.

The measuring ranges refer to the system and are defined by the technical limits of the meter and test strip. The measuring ranges are as follows: INR: 0.8-8.0 / % Quick: 120-5 / Seconds: 9.6-96

The therapeutic range is defined by your treating physician, who will tell you what your individual therapeutic range is. Your test result should fall within the therapeutic range.

⚠ **How to proceed if test results are outside the measuring range or therapeutic range**

If the meter displays < (less than) 0.8 or > (greater than) 8.0 INR, repeat the test. If, after repeating the test, you get the same display (either < 0.8 or > 8.0), call your doctor.

If the measured result is outside the therapeutic range specified by your treating physician, repeat the test. If the result is still outside the therapeutic range immediately contact your physician and ask for the appropriate measures to take.

Notes on PT self-testing

PT self-testing supplements physician care, but is not a substitute for it. PT self-testing gives the anticoagulated patient more security in everyday life. Results should be recorded in a notebook (patient’s logbook) and presented to the physician at every consultation. This makes it easier for the physician to assess the overall quality of therapy management.

Storage and stability

Store at 2-30 °C.

The test strips can be used up until the expiry date printed on the box and test strip container.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

⚠ This is required in order to prevent the remaining test strips from deteriorating through exposure to external influences like humidity.

Materials provided

- Test strips and 1 code chip

Materials required (but not provided)

- [REF] **07404379**, CoaguChek INRange meter or
- [REF] **04625412**, CoaguChek XS meter
- Lancing device, e.g. CoaguChek XS Softclix
- Lancoets, e.g. CoaguChek Softclix lancet

⚠ The CoaguChek XS Softclix lancing device has been specifically developed for self-testing. It is not suitable for testing different persons as there is a risk of contamination.

CoaguChek INRange

CoaguChek XS

Sample material

Use only fresh capillary blood.

Before you test

Before carrying out your first blood test, please read the meter manual for your CoaguChek INRange / CoaguChek XS meter carefully to acquaint yourself with the way it works.

Getting a good drop of capillary blood

Increasing the blood flow in your finger will help you to get a good drop of blood. Before you lance your finger, try the following techniques until you see that your fingertip has good color.

- Warm your hand, e.g. by washing your hand with warm water.
- Hold your arm down to the side so that your hand is below your waist.
- Massage your finger from its base.

Performing the test

⚠ Make sure that the code chip number matches the code number printed on the test strip container. If you use a new test strip lot you first have to insert the corresponding code chip provided with the strips, otherwise you cannot perform a test.

- Wash your hands with warm, soapy water. Dry completely. Residues of water on the skin can dilute the drop of blood and produce false results.
- Place the meter on a level, vibration-free surface or hold it in your hand so it is roughly horizontal.
- Remove 1 test strip from the container and immediately close the container. Make sure it seals tightly. Do not touch the test strip with wet hands. This may damage the test strips. Use the test strip within 10 minutes after taking it out of the container.
- Hold the test strip so the arrows and blood drop symbols are facing upward. Insert the test strip as far as you can into the test strip guide in the direction indicated by the arrows. The meter powers ON.

- CoaguChek INRange:** The meter reads the strip information and checks if it matches the code chip already stored in the meter. Only in case of a match, the meter will automatically proceed.

CoaguChek XS: The meter reads the code chip information and displays the code number on the screen. If the number displayed matches the code number printed on your test strip container, press the M button to proceed.

CoaguChek INRange	CoaguChek XS
code match: meter proceeds	code match: press M button to proceed

In case of a mismatch please insert the correct code chip provided with the test strips.

- The meter warms up, which takes up to 30 seconds.
- When the meter is warmed up, a flashing test strip and blood drop symbol appear and the meter begins a countdown. You have 180 seconds to apply blood to the test strip.
- Use the lancet device to prick your finger. Gently massage the lanced finger until a drop of blood is formed. Do not press or squeeze the finger. See the meter manual for more information.
- Apply the blood drop to the test strip within 15 seconds after lancing your finger. Applying blood after this period of time would falsify the result. If it takes longer to form a good drop of blood, lance a different finger for doing another test.
- Apply the blood to the side of the test strip or from the top. Only the blood drop should make contact with the strip, avoid touching the strip with the skin. It's important to hold the blood drop to the test strip until you hear a beep.
- Then, pull your finger away from the test strip. Do not touch the test strip until the result is displayed.
- The result appears within 1 minute.

ⓘ **Make sure you are using the code chip that came with the test strips. Use the test strip within 10 minutes from the time you took it out of its container. Apply the blood drop to the test strip within 15 seconds after lancing. Remember to apply only one drop of blood—don't add more. Do not touch or remove test strip when a test is in progress.** Record the test results in your patient logbook and present this to your physician at every consultation. Results are also saved automatically to the meter’s memory. After testing, you can discard the used lancet and test strip with your regular household waste. Dispose of used lancets carefully (e.g. use a sturdy sharps container with lid). If you need to redo a test, use a new lancet, a new test strip, and a different finger.

Quality control

The CoaguChek INRange and the CoaguChek XS system have a number of inbuilt quality control functions. The meter performs the quality control testing as part of every blood test. For more information on the built-in quality control functions, please refer to the according meter manual.

Cleaning and disinfection instructions

To avoid malfunction of the meter, please follow the instructions provided in the according meter manual section.

Possible causes of errors

The CoaguChek INRange and the CoaguChek XS meter continually check their systems for unexpected and unwanted conditions. An error message appears if the meter has a problem that needs your attention.

For CoaguChek XS: If you receive an error message, please refer either to the meter manual or contact your local customer support and service center.

For CoaguChek INRange: Read the error message and take the action suggested in the message.

If problems occur during testing, please check the following:

- The blood drop must be a minimum of 8 µL in volume. Low sample volume will cause an error message.
- Are you testing exactly as described in the according meter manual and also in the "Getting a good drop of capillary blood" section in this document? Please read the instructions carefully.
- Always make sure the test strips are correctly stored (see section "Storage and stability") and that you carry out the test within 10 minutes of removing a test strip from its container.

- Please check date and time** on the meter display. Incorrect settings might cause an error message.
- Is the test strip guide dirty? Clean the meter as described in the meter manual. Repeat the test with a new test strip.

ⓘ **If the problem persists, please contact your local Roche Diagnostics customer support (for contact details see your meter manual).**

⚠ Certain errors might occur sporadically due to an activation of the system fail safe mechanisms that are designed to prevent the release of wrong measurement results. In rare cases, these errors might also be received by **patients under clinical conditions**, e.g. treatment with vitamin K antagonists in combination with antibiotics and/or chemotherapeutics, or if the blood sample contains extremely high concentrations of oxidizing substances, e.g. after a vitamin C infusion. Also, in rare cases, **patients with abnormal or unusually long clotting times** may receive these error messages.

CoaguChek INRange	CoaguChek XS
Error "M-44"	"error 6"
Error "W-45"	"error 7"

ⓘ **Repeat test with a new strip. If the error is displayed repeatedly, you must use an alternative test method to confirm the result. Please contact your physician without delay.**

Additional information for Health Care Professionals

Reagents

The test strip contains human recombinant tissue factor as activator, stabilizers and preservatives.

Precautions and warnings

For in vitro diagnostic use. Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Limitations - interference

Testing performed with the following in-vitro spiked samples or native blood samples indicated no significant effect on test results:

- Ascorbic acid up to 170 µmol/L (30 mg/L)
- Bilirubin up to 513 µmol/L (30 mg/dL)
- Hemolysis up to 0.62 mmol/L (1000 mg/dL)
- Triglycerides up to 5.7 mmol/L (500 mg/dL)

Hematocrit ranges between 25 % and 55 % do not significantly affect test results.

The CoaguChek XS PT Test PST is insensitive to unfractionated and fractionated heparin concentrations up to 2 IU/mL blood.

Information about pharmacokinetics and the dosage-plasma-relation of the respective heparin are provided in the package insert by the legal manufacturer of the drug.

PT results may be influenced from patients treated with Direct Oral Anticoagulants (DOACs) e.g. rivaroxaban, apixaban, edoxaban, betrixaban and dabigatran or anticoagulants other than vitamin K antagonists (e.g. hirudin and other thrombin inhibitors). For such patients, medical decisions should not be based on CoaguChek system measurements.

Note: Samples from patients treated with the following drugs must not be tested with the system: protamine sulfate, oritavancin, calcium dobesilate, fondaparinux.

⚠ **The action of oral anticoagulants (coumarin derivatives) can be increased or weakened when other medication is taken simultaneously** (e.g. antibiotics, but also prescription-free medication like pain relievers, antirheumatic medication and medication against influenza). This, in turn, can also lead to either an increase or a decrease in prothrombin time (INR). If other medication is taken, it is recommended that the prothrombin time be checked more frequently and that the anticoagulant dose be subsequently adjusted.

Antibodies against domain I of β2-glycoprotein I, also known as Anti-phospholipid antibodies (APA) or Lupus anticoagulant (LA), may prolong PT.

Unexpected PT values should always be followed by additional testing to determine the source of the influence.³

Drug interferences are measured based on recommendations given in CLSI guidelines EP07 and EP37 and other published literature. Effects of concentrations exceeding these recommendations have not been characterized.

For diagnostic purposes, the results should always be assessed in conjunction with the patient’s medical history, clinical examination and other findings.

Specific Performance Data

Equivalence of CoaguChek XS and CoaguChek INRange was shown. Data presented here were obtained with the CoaguChek INRange system and are representative for both systems.

Precision

Repeatability of CoaguChek XS PT Test PST was determined with duplicates of fresh capillary whole blood samples at 4 external sites.⁴

Repeatability			
Range (INR)	Number of duplicates	SD (INR)	CV (%)
< 2.0	194	0.04	3.8
2.0-3.0	326	0.07	2.9
3.1-4.5	138	0.12	3.4
> 4.5	28	0.08	1.5

Method comparison

INR based

A method comparison at 3 external sites was performed with CoaguChek XS PT Test PST using fresh capillary whole blood versus the reference method Innovin

(Siemens) using venous citrate blood collected with Sarstedt blood collection tubes.⁵

Number of samples measured: 1020

Passing/Bablok⁶

y = 1.07x - 0.04 INR

Kendall's τ = 0.86

The PT values were between 0.9 and 5.2 INR.

Note: While the INR unit was created in general for the best possible comparison between methods (Point-of-Care vs. laboratory methods and between laboratory methods), several factors can have a significant or even systematic influence on the comparison of PT/INR results obtained with different methods. The most important factor, when changing methods, is the type of thromboplastin used (i.e. human recombinant, rabbit or bovine).⁷

The CoaguChek method uses human recombinant thromboplastin. Therefore, the comparability to tests using other human recombinant thromboplastins is best, whereas higher deviations can occur with other thromboplastin types.

However, those higher differences between thromboplastins of different (rabbit, bovine) origin are not an issue specific for CoaguChek assays. Similar differences can be observed when a human recombinant thromboplastin-based laboratory method is compared with several other (rabbit, bovine) laboratory methods. To minimize these differences, in a monitoring situation, it is recommended that each laboratory uses results from a method using only one type of thromboplastin for each patient.

If testing methods for the comparison of patient INR values are changed, especially when methods with thromboplastins of a different origin are used, the deviations which may occur need to be taken into account.

Other potential influencing factors:

- Variations between different laboratory instruments and between different reagent lots
- Pre-analytics, i.e. sample tubes of different manufacturers used for laboratory testing

Sample collection tubes other than Sarstedt for the PT INR determination with the reference method may lead to different results in method comparison.^{8,9}

For further information, please refer to the appropriate meter manual and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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Letteratura / Literatur

- Leichsenring I, Plesch W, Unkrig V, et al. Multicentre ISI assignment and calibration of the INR measuring range of a new point-of-care system designed for home monitoring of oral anticoagulation therapy. Thromb Haemost 97:846-861, 2007.
- Plesch W, Van den Besselaar AMHP: Validation of the international normalized ratio (INR) in a new point-of-care system designed for home monitoring of oral anticoagulation therapy. Int Jnl Lab Hem 31:20-25, 2009.
- Luxembourg B, Plesch W, Lindhoff-Last E. INR results of the CoaguChek® XS system in patients with lupus anticoagulants agree with Innovin INR. Poster presented at the 52nd Annual Meeting of the GTH (Gesellschaft für Thrombose- und Hämostaseforschung) in Wiesbaden, Germany, Feb 20–23, 2008.
- Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2014.
- Multicenter Performance Evaluation of the CoaguChek INRange System; White Paper from Roche Diagnostics available upon request.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem. 26(11):783-90,1988.
- Plesch W, Wolf T, Breitenbeck N et al. Results of the performance verification of the CoaguChek XS System. Thromb Res 2008;123:381-389.
- van den Besselaar AMPH, Hoekstra MMCL, Witteveen E, et al. Influence of Blood Collection Systems on the Prothrombin Time and International Sensitivity Index Determined With Human and Rabbit Thromboplastin Reagents. Am J Clin Pathol 127:724-729, 2007.

- van den Besselaar AMPH, Rutten WPF, Witteveen E. Effect of magnesium contamination in evacuated blood collection tubes on the prothrombin time test and ISI calibration using recombinant human thromboplastin and different types of coagulometer. Thrombosis Research 115:239-244, 2005.

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© 2020, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim www.roche.com